



SMEARS

H.P. Techs prefer our cotton smear !!! Why ???

Uniformity – each smear is the same. They are made of 8 ½ oz 100% cotton twill. They do not curl when used on a rough surface. Also, the adhesive leaves no residue in the planchet.

We have modified this smear to help **reduce** it's cost. Our record folder has been down sized to 2 5/8" W x 3 3/4" H. For your convenience the smear and record folder are packed loose and unfolded. The box is plain white, but sturdy enough for stacking and rough handling.

Please note that the quality of the smear itself has not been altered.

DATE _____	TIME _____
LOCATION _____	
_____ NO. _____	
TECHNICIAN _____	
COUNT	TYPE
..... FOLD HERE	
800-933-8501	G/O Corporation

Ordering Information:

GJ1130-250 – Cotton Twill with record folder – 250 per box

GJ1130-500 – Cotton Twill with record folder – 500 per box

GJ1075 Disc 500 – Cotton Twill – disc only – no record paper – 500 per box

GJ1302 500/box - Cotton Twill - Piggy Back





TECHNICAL DATA

	IMPERIAL	METRIC
Adhesive Thickness:	1.0 mils (0.0010")	0.0254mm
Peel Adhesion (PCTC-1)	16 oz/inch width	4.5 N/2.5 cm
Shear Adhesion (PSTC-7)	12 hours @ 1.1 psi	12 hours @ 7.6 kPa
Maximum Continuous Use Temperature:	225 degrees F	107 degrees C
Minimum Application Temperature:	0 degrees F	-18 degrees C
Release Liner:	70lbs / ream Polycoated bleached kraft	115 GSM

The information listed is in good faith and has been established from sources believed to be secure and reliable. The values listed are typical properties and are not intended to be used as specification for defined applications. User must determine the product suitability for the end us application and assume all risks and liabilities.

